

KSEP20L-CAN

Cesana AG

Baiergasse 56, CH-4126 Bettingen, Switzerland

www.cesana-ag.ch info@cesana-ag.ch

REMBE® Research+Technology Center GmbH

Zur Heide 39, D-59929 Brilon, Germany

www.rembe-rtc.de info@rembe-rtc.de

History



The 20-liter apparatus was developed in collaboration with Ciba-Geigy (Lütolf, Bartknecht, Siwek). This is why it is also known today as the Siwek sphere. The test method was developed by Ciba-Geigy and the apparatus by Kühner AG (Cesana). Our 20-liter apparatus quickly found its way into the following standards:

ASTM E1226: Explosibility of Dust Clouds

ASTM E1515: Minimum Explosible Concentration MEC of Combustible Dusts

ASTM E2931: Limiting Oxygen Concentration LOC of Combustible Dusts

EN 14034-1: Determination of max. explosion pressure P_{max} ...

EN 14034-2: Determination of max. rate of explosion pressure rise $(dP/dt)_{max}$...

EN 14034-3: Determination of the lower explosion limit LEL ...


EN 14034-4: Determination of the limiting oxygen concentration LOC ...

EN ISO/IEC 80079-20-2: Material characteristics. Combustible dust methods

Kühner AG discontinued the production of safety testing equipment at the end of 2016. These products will be continued by Cesana AG and, from 2025, together with REMBE® Research+Technology Center GmbH.


Calibration-Round-Robin CaRo

According to the standards, the 20-liter apparatus must be calibrated at regular intervals. For this reason, we have been carrying out an international calibration round robin (CaRo) periodically since 1993. Our experience with the calibration round robin was the guideline for the development of the new 6th generation instrumentation: improved test quality and simplified operation.



New instrumentation for 20-liter apparatus

Over the almost 5 decades, the instrumentation has always been adapted to the state of the art. The aim has always been full backwards compatibility. Equipment from 1976 onwards can be brought up to date.



Principle of new instrumentation

Absolute pressure measurement


The explosion parameters P_{max} and K_{max} are proportionally influenced by the initial pressure P_i , the pressure at the time of ignition. According to the standards: $P_i = 1013 \text{ mbar} = \text{normal pressure}$.

The scattering of P_{max} in the image "Pmax of CaRo23" is mainly due to deviations from P_i .

Instead of the analog and relative pressure display for sphere and dust storage, **absolute** and **digital** measurements are now used. Errors caused by the barometric ambient pressure are eliminated.

Evacuation of the sphere and filling the dust storage take place automatically

The results of the preparation are displayed in addition to the explosion characteristics. The absolute initial pressure P_i is calculated from the measured values of P_v and P_z .



| dP/dt | | | |
|--------------------------|----------|-----------|--------|
| Test: 49 (3) 125 [g/m3] | | | |
| | setpoint | effective | FS [%] |
| P_v vacuum [mbar] | 413 | 410 | -0.3 |
| P_z air [bar] | 21.01 | 21.07 | 0.3 |
| P_i calculated* [mbar] | 1013 | 1012* | |


Dynamic pressure measurement

Although the pressure curves are still measured relatively using the robust and proven Kistler piezo pressure transducers, the absolute initial pressure is known and is taken into account during recording. The new Kistler charge amplifiers have a digital design and communicate with the KSEP333 via IO-Link.

As before, the P_{ex} and P_m values are given as gauge pressure. However, they can be related to a normal pressure of 1013 mbar.

Pressure rise dP/dt

For the decision "ignition / no ignition" for LEL and LOC, the representation of dP/dt can be useful (see ASTM standards) and helps to differentiate between "igniter" and "dust explosion".



Limiting Oxygen Concentration (LOC)


Thanks to the precise and absolute pressure measurement in the sphere and dust storage, it is possible to produce the desired air/inert gas mixture using the partial pressure method.

Only one inert gas cylinder is required. Favorable inert gas cylinders contain a residual oxygen content. This is taken into account in the calculation.

The gas mixtures are created automatically.


Example:

Set point $O_2 = 9\%$ (Inertgas oxygen = 1%)




Sphere (Pv)

1. evacuate




2. fill inertgas




Dust storage (Pz)

3. fill air



4. fill inertgas



Inaccuracies during preparation and leaks in the sphere can influence the mixture. The effective gas mixture can be calculated by precisely measuring the absolute pressures in the sphere and dust storage.

In contrast to the previous partial pressure method, the same mixture is now prepared in the dust storage and sphere.

This also applies to the "hybrid mixture" and "gas" test methods.

| dP/dt Test: 12 (2) 60 [g/m3] | | | |
|------------------------------|----------|-----------|--------|
| | setpoint | effective | FS [%] |
| Pv Vakuum [mbar] | 166 | 168 | 0.2 |
| Pv Inertgas [mbar] | 413 | 412 | -0.1 |
| Pz Luft [bar] | 21.01 | 21.02 | 0.0 |
| Pz Inertgas [bar] | 13.66 | 13.66 | 0.0 |
| Pi berechnet* [mbar] | 1013 | 1012* | |
| O2 berechnet* [vol%] | 9.0 | 9.0* | |
| Inertgas O2 [vol%] | | 1.0 | |

Temperature control

According to ISO/IEC 80079: Initial temperature = $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$

The jacket temperature is continuously measured, controlled and transferred to the table for each test.


Connection for an external solenoid valve for simple on/off control of the cooling water flow and thus the jacket temperature. This significantly reduces the consumption of cooling water.

Safety

By monitoring the position of the two ball valves (outlet and vacuum), the operator is guided through the test sequence and protected against incorrect operation.

Calibration of static pressure measurement „Pv“

When installing the software, the altitude above sea level of your device must be entered. The current static pressure is then converted to sea level and displayed. Calibration is carried out by comparison with the current meteorological data.



Calibration of pre-chamber pressure „Pz“


After an empty test without dust and igniter, $P_i = 1013 \text{ mbar}$ = normal pressure must be present in the sphere. This is measured with the static pressure measurement (Pv). Calibration chain: Meteodata → Pv → blank test → Pz

Calibration of dynamic pressure measurement

The sensitivity of the pressure sensors [pC/bar] is set digitally via IO-Link.

The system consists of 2 independent measuring channels. As long as both match, it can be assumed that the measured values are correct. If they do not match, the charge amplifier can be easily checked:

A generator for testing with a rectangular charge signal is built in.




New control parameter „tdi“

Time delay of the ignitors = time difference between the electrical energization of the ignitors and the first pressure raise.

Carbonized coating on the electrode rods is a poor conductor and reduces the current flow to the ignitors. As a result, the ignition pills are activated with a delay. If the maximum permissible time delay is exceeded, an error message is issued.

Unfortunately, this error is often observed when evaluating CaRo ring tests.



| Technical data: | 20-I-Sphere | Control unit KSEP333 | |
|-------------------------------|---------------|-----------------------------|-------------------|
| Material of construction No.: | 1.4435 | PC-Software | KSEP 8.0 |
| Volume of sphere: | 20 liters | Dispersion pressure (Pz) | absolute, digital |
| Volume of water jacket: | 1,5 liters | Pressure in the sphere (Pv) | absolute, digital |
| Design pressure of sphere: | 30 bar (60°C) | Measurement pressure curve | relative, digital |
| Design pressure of jacket: | 10 bar (60°C) | Temperature measurement | digital |
| Test pressure of sphere: | 42.9 bar | Charge amplifiers | Kistler 5028A |
| Test pressure of jacket: | 15.8 bar | Interface to PC | CAN - USB |
| Design temperature: | 60 °C | Mains connection | 100 ... 240VAC |
| Bayonet ring aperture Ø: | 96 mm | | |
| Cleaning aperture Ø: | 140 mm | | |
| Sight glass Ø: | 30 mm | | |

Order data

| | |
|---------|---|
| SP3003 | Complete test equipment KSEP20L-CAN |
| SP30031 | Upgrade: KSEP20L → KSEP20L-CAN comprising: Electronics unit KSEP333-CAN Pneumatic unit KSEP311 Accessories |
| SP30333 | Upgrade: KSEP33x → KSEP333-CAN Without pneumatic unit KSEP311. (No absolute pressure measurement) |

References

| | |
|---|---|
| Bartknecht W. Zwahlen G. | 1993, Explosionsschutz: Grundlagen und Anwendung Springer, Berlin, DE |
| Siwek. R. | Explosion Characteristics and Influencing Factors, International Symposium on Control & Prevention of Dust Explosions, Basel, 1982 |
| Siwek. R. | Reliable Determination of the Safety Characteristics in 20-I-apparatus. Symposium on Flammable Dust Explosion; St. Louis Missouri, USA, 1988 |
| Cesana C. Eiche M. Schwaninger M. | Quality Management in the Determination of Safety Characteristics Chemical Engineering Vol.75, 2019 |
| Cesana C. Siwek R. | Manual for the 20-I-apparatus 8.0, 2024 |